

# Eurotrode

## ELECTRODES SPECIFIC FOR ANY NEED!

 **Pirrone** S.r.l.  
ELECTRODES

### FOR SHORT, MEDIUM TERM MONITORING AND RESTING PROCEDURES



#### PLG50

**Size:** Ø50mm  
**Back:** foam  
**Gel:** liquid  
**Sensor:** Ag/AgCl  
60 pcs per pouch  
600 pcs per box  
7200 pcs per carton

Re-order code: 12875



#### PSG30

**Size:** 35mm x 27mm  
**Back:** foam  
**Gel:** solid  
**Sensor:** Ag/AgCl  
100 pcs per pouch  
1200 pcs per box  
14400 pcs per carton

Re-order code: 12777



#### PLG30 (pediatric)

**Size:** 35mm x 27mm  
**Back:** foam  
**Gel:** liquid  
**Sensor:** Ag/AgCl  
60 pcs per pouch  
600 pcs per box  
7200 pcs per carton

Re-order code: 12909



#### PSGST-064M

**Size:** 44mm x 28mm  
**Back:** foam  
**Gel:** solid  
**4mm banana connector**  
60 pcs per pouch  
600 pcs per box  
7200 pcs per carton

Re-order code: 12773



#### PLG45

**Size:** 45mm x 42mm  
**Back:** foam  
**Gel:** liquid  
**Sensor:** Ag/AgCl  
60 pcs per pouch  
600 pcs per box  
7200 pcs per carton

Re-order code: 12987



#### PFR2034 TAB ELECTRODE

**Size:** 34mm x 20mm  
**Back:** foam  
**Sensor:** Ag/AgCl  
100 pcs per pouch  
5 pouches per box

Re-order code: 12774

Easy to apply and remove - Highly conductive and hypoallergenic gel - Latex free - Not irritating - Strong but skin-friendly adhesive  
Meet and exceed all ANSI/AAMI performance specifications

## FOR STRESS TEST, HOLTER AND LONG-TERM MONITORING



Re-order code: 12874

### PSG50

**Size:** Ø50  
**Back:** foam  
**Gel:** solid  
**Sensor:** Ag/AgCl  
60 pcs per pouch  
600 pcs per box  
7200 pcs per carton



Re-order code: 12988

### PSG50C

**Size:** Ø50  
**Back:** TNT  
**Gel:** solid  
**Sensor:** Ag/AgCl  
60 pcs per pouch  
600 pcs per box  
7200 pcs per carton



Re-order code: 12979

### PSG51

**Size:** 51mm x 33mm  
**Back:** foam  
**Gel:** solid  
**Sensor:** Ag/AgCl  
100 pcs per pouch  
1000 pcs per box  
12000 pcs per carton



Re-order code: 12958

### PSG45

**Size:** 45mm x 42mm  
**Back:** foam  
**Gel:** solid  
**Sensor:** Ag/AgCl  
60/100 pcs per pouch  
600/1000 pcs per box  
7200/12000 pcs per carton



Re-order code: 12876

### P40SG

**Size:** 40mm x 36mm  
**Back:** foam  
**Gel:** solid  
**Sensor:** Ag/AgCl  
60/100 pcs per pouch  
600/1200 pcs per box  
7200/14400 pcs per carton

## LOW-COST AND MULTIFUNCTION ELECTRODES



Re-order code: 12978

### PLG44ST

**Size:** 44mm x 32mm  
**Back:** foam  
**Gel:** liquid  
**Sensor:** Ag/AgCl  
60 pcs per pouch  
600 pcs per box  
7200 pcs per carton



Re-order code: 12973

### PSG44ST

**Size:** 44mm x 32mm  
**Back:** foam  
**Gel:** solid  
**Sensor:** Ag/AgCl  
100 pcs per pouch  
1200 pcs per box  
14400 pcs per carton

Easy to apply and remove - Highly conductive and hypoallergenic gel - Latex free - Not irritating - Strong but skin-friendly adhesive  
Meet and exceed all ANSI/AAMI performance specifications

## RADIOTRANSSPARENT FOR RMN AND X-RAY PROCEDURES



Re-order code: 12986

### PLG45RT

**Size:** 45mm x 42mm  
**Back:** foam  
**Gel:** liquid  
**Sensor:** carbon  
60 pcs per pouch  
600 pcs per box  
7200 pcs per carton



Re-order code: 12985

### PSG45RT

**Size:** 45mm x 42mm  
**Back:** foam  
**Gel:** solid  
**Sensor:** carbon  
100 pcs per pouch  
1000 pcs per box  
12000 pcs per carton




Re-order code: 12984

### PSG50SRT

**Size:** Ø50mm  
**Back:** foam  
**Gel:** solid  
**Sensor:** carbon  
60 pcs per pouch  
600 pcs per box  
7200 pcs per carton


## BBYTRODE NEONATAL ELECTRODES



Re-order code: 12995

### PN378

**Size:** 22mm x 22mm  
**Back:** TNT  
**Gel:** solid  
**Connector:** 1,5mm DIN  
3 pcs (blister) per pouch  
100 pouches per box  
12 boxes per carton




Re-order code: 12996

### PN378RT

**Size:** 22mm x 22mm  
**Back:** TNT  
**Gel:** solid  
**Connector:** 1,5mm DIN  
(with radiotrasparent wires for RMN)  
3 pcs (blister) per pouch  
100 pouches per box  
12 boxes per carton


# EUROCLAMPS AND EUROSUCTIONCHEST FOR SHORT MONITORING AND RESTING PROCEDURES



## Euroclamp for adult

Sensor: Ag/AgCl  
4pcs per pouch

Re-order code: 12917



## Euroclamp pediatric

Sensor: Ag/AgCl  
4pcs per pouch

Re-order code: 12918



## Eurosuctionchest 15mm

Ag/AgCl ring 15mm  
4mm socket and snap  
6 pcs per pouch

Re-order code: 12919



## Eurosuctionchest 24mm

Ag/AgCl ring 24mm  
4mm socket and snap  
6 pcs per pouch

Re-order code: 12920



## Eurosuctionchest 30mm

Ag/AgCl ring 30mm  
4mm socket and snap  
6 pcs per pouch

Re-order code: 12922

CERTIFICATES IN FORCE



Certificate n.  
9124.PIR2  
EN ISO  
13485:2012



Certificate n.  
9120.PIRR  
ISO  
9001:2008



REV. 01/2018

Universal connector for any connection type

 **Pirrone** S.r.l.  
ELECTRODES

Pirrone S.r.l. • via Balduccio da Pisa n. 12 • 20139 Milano • Italy  
tel. [+39] 02.525397.1 • fax [+39] 02.57404077  
<http://www.pirrone.com> • e-mail: [info@pirrone.com](mailto:info@pirrone.com)